



SINGLE USE ELASTIC STAY RETRACTORS FOR ENT

SPECIFICALLY DESIGNED AND PACKAGED FOR HEAD AND NECK PROCEDURES

Bailey

INSTRUMENTS LTD

Bailey Instruments Ltd

Unit 19, Guinness Road Trading Estate,
Guinness Rd, Trafford Park, M17 1SB

Tel: 0161 872 8707 | Email: sales@baileyinstruments.co.uk
www.baileyinstruments.co.uk



MD642415



441219

Bailey

INSTRUMENTS LTD

Aseptically Packaged 4 stay pack now available via the NHS Catalogue / NHS Supply Chain

The system delivers precise location and retraction force in soft tissue. This can create a stable operative field with minimal field movement during surgery.

Contact us for further information.

Clinical Applications:


Skin Incision/Closure, Tracheostomy, Lymph Node Biopsy, Excision of Thyroglossal Duct Cyst, Parotidectomy, Thyroidectomy and Parathyroidectomy, Lateral Neck Access. Designed for use with Retractor Rings, Anchor Plates or Artery Forceps.

Technical Specification and ordering details overleaf.

TECHNICAL SPECIFICATION – SUSOL® ELASTIC RETRACTOR STAYS

MPC	NPC	SECONDARY DESCRIPTION	SALES UNITS	NHSSC SELL PRICE (INC VAT + DEL)
SS25-6001	FFQ4213	susol® Retractor Elastic Stay - Hook 5mm	16 stays - 4 x 4 stay packs	£151.30

- **Materials:** Stainless Steel Hook embedded in silicone.
- **Performance:** Stress tested for retraction tension and stability
- **Packaging:** Aseptic presentation in easy access plastic tray
- Material 0.5 PET
- Manufactured under the requirements of the Labour Standards Assurance System
- Latex Free



More product developments to follow for Plastics, ENT and Max Fax. Give us your ideas and we'll turn them into products.

